



THOMAS MEDICAL

SHAPEABLE HS CATHETER



SHAPEABLE HS CATHETER

- Sterile, Disposable, Ready to use
- 5Fr and 7Fr sizes available
- Integrated stylet is malleable and can be re-shaped as necessary
- Angle can be set to align with cervical canal
- Stiffness beneficial for restricted or stenotic cervical canal
- For use with saline or water-based contrast media
- THIS PRODUCT DOES NOT CONTAIN LATEX

TMI1185 5Fr HS Catheter, Box of 10

TMI1187 7Fr HS Catheter, Box of 10

THOMAS MEDICAL

5610 W 82ND STREET

INDIANAPOLIS, IN 46278

+1.800.556.0349 TF

+1.888.411.3754 FAX

WWW.THOMASMEDICAL.COM