



THOMAS MEDICAL

NOW
FEATURING
DIRECTIONAL
INJECTION

IUI CATHETER

A Catheter For Intrauterine Insemination



IUI Catheter



Shapeable IUI Catheter

IUI Catheter

- High quality IUI catheters designed for intrauterine insemination
- Atraumatic tip designed to facilitate smooth insertion
- Shapeable IUI catheter has built-in stylet which is malleable and can be shaped as necessary
- Depth marks are indicated on both catheters for accurate placement
- Side hole at tip is aligned with the printed depth marks. Clinician can align catheter to inject toward the left or right fallopian tube ostium. This feature is shown below:



TMI1192 - IUI CATHETER, 15 CM

TMI1193 - SHAPEABLE IUI CATHETER, 15 CM

TMI1194 - IUI CATHETER, 18 CM

TMI1195 - SHAPEABLE IUI CATHETER, 18 CM

TMI1196 - IUI CATHETER, 23 CM

TMI1197 - SHAPEABLE IUI CATHETER, 23 CM

THOMAS MEDICAL

5610 W. 82ND STREET
INDIANAPOLIS, IN 46278

(800) 556-0349 PH
(888) 411-3754 FX

WWW.THOMASMEDICAL.COM

IDEAL FOR: • Intrauterine Insemination